

QUICKSTART GUIDE

RESCUE[®] Concentrate


Rescue[™] Concentrate Dilution and Contact Times Overview

DAILY USE	BACTERICIDAL FUNGICIDAL (PARVO)	VIRUCIDAL	NON-FOOD CONTACT SURFACE SANITIZER
5 Min @ 1:64 2 oz/gal	5 Min @ 1:16 8 oz/gal	5 Min @ 1:64 2 oz/gal	3 Min @ 1:128 1 oz/gal
Surface should be visibly wet for contact time.			

Rescue[™] Concentrate can be diluted at specific ratios based on the type of application you require. Excellent for daily use and outbreak situations. Rescue[™] Concentrate has fast contact times. Make sure that the surface is visibly wet for the entire contact time.

For more information on a specific organism, use the Pathogen Claim Tool on the Rescue[™] Concentrate page at www.LearnAboutRescue.com.

How to Use Rescue[™] Concentrate



STEP 1 Dilute the concentrate


Use a mixing system, approved accessories, or a measuring cup. For every day use, dilute the concentrate at a 1:64 ratio (2 oz per gallon) but be sure to follow the proper dilution measurements for specific pathogens.

Approved Accessories: Watch Video on www.LearnAboutRescue.com



STEP 2 Remember to apply a workplace label to the solution container

Use a Rescue[™] Concentrate workplace label and record the dilution ratio, the lot number and expiry date.




STEP 3 Pre-clean heavily soiled areas

Remove visible soil and pre-clean the area by working from the least contaminated area to the most contaminated area.

STEP 4 Apply the diluted solution until thoroughly wet

The entire surface area should be covered and appear visibly wet. Wettability and surface coverage are equally as important to ensure disinfection can be achieved in one application.



STEP 5 Allow surface to remain visibly wet for 5 minutes

In order for disinfectants to be effective they must remain wet on the surface for their entire contact time. Be sure to follow the proper contact times for specific pathogens.

STEP 6 Wipe surface dry or allow to air dry

Use a damp, clean cloth to dry small areas, or for larger areas like floors let them air dry and rinse if needed.

Success & Troubleshooting Tips

Start Fresh.

If you are switching from a previous cleaner and disinfectant to using Rescue[™] for the first time, start with fresh accessories and tools. Rescue[™] will remove build up from previous disinfectants and may cause a non-toxic reaction and musty odor. The odor will eventually disappear after a few times of using Rescue[™].

Measure. Measure. Measure.

There are many mixing systems available, but at the very minimum use a measuring cup! Measuring will ensure proper dilution and ensure proper claims are being met when you want to kill certain pathogens. Use test strips to confirm that you have mixed the correct dilution.

More is NOT Better.

Follow the dilution ratios and you will conserve product and prevent unnecessary damage and sticky floors. When diluted properly, Rescue[™] Concentrate is a cost-effective solution for your everyday cleaning and disinfectant needs.

Wait Before Wiping Dry.

Always follow the required contact time when disinfecting a surface for a particular pathogen. Consider setting a timer and walking away until the timer goes off.

For more information visit:
Rescuedisinfectants.com

Virox Animal Health
2770 Coventry Road
Oakville, ON L6H 6R1 CANADA
Toll Free: 1-800-387-7578